

Infrared Analysis of Thin Multilayered Polymer Film Using Cantilever Enhanced Photoacoustic Detector

Kaori Watanabe and Fumiko Hiwatashi  Systems Engineering

Systems Engineering Inc. www.systems-eng.co.jp info@systems-eng.co.jp

Introduction

Polymer thin film is used in...

- Structural material
- Electronic device
- Adhesive
- Paint and printing
- Biomaterial

nm to μm ordered polymer structure
Multilayered films for hybrid and novel functions
Depth profiling analysis is necessary


Surface analysis

- XPS, SSIMS, SPM, XRD etc.
- Vibrational spectroscopy -> definitive molecular information
» PAS: Non-destructive (little sample preparation), depth-profiling

Novel PAS detector (GASERA; booth 555)



- Cantilever and laser interferometer based detector
- **Ten times more sensitive than traditional acoustic microphone**



This study

- Infrared PAS analysis of multilayered films
- Determination of accurate sampling depth of PAS with cantilever detector

Photoacoustic Spectroscopy (PAS)


Photoacoustic spectroscopy

1. Sample absorbs incident light
2. Absorbed energy partially transformed into heat (non-radiating relaxation)
3. Heat diffuses back to sample surface
4. Thermal expansion of surrounding gas detected

PAS advantage

- Highly sensitive
- Non-destructive
- Depth-profiling
- Versatile
 - » Irregular shaped, small, or opaque sample etc.
 - » THz, IR, NIR, UV-Vis etc.

PAS Sampling Depth




$$\begin{aligned} \alpha: & \text{Thermal diffusivity } \alpha = k/pC [\text{cm}^2 \text{s}^{-1}] & \omega: & \text{Angular modulation frequency } \omega = 2\pi\nu \\ k: & \text{Thermal conductivity } [\text{W cm}^{-1} \text{K}^{-1}] & V: & \text{Mirror velocity } [\text{cm s}^{-1}] \\ p: & \text{Density } [\text{kg cm}^{-3}] & v: & \text{Infrared wavenumber } [\text{cm}^{-1}] \\ C: & \text{Specific heat } [J \text{ kg}^{-1} \text{ K}^{-1}] \end{aligned}$$

Results and Discussion

1. PAS spectrum of multilayered film

Layered sheets of food wrap
PE 10 μm , PVC 10 μm , PVDC 10 μm
Mirror velocity $V = 5 \text{ kHz} = 0.32 \text{ cm s}^{-1}$




1740 cm^{-1} peak of PVC observed

Sampling depth $> 10 \mu\text{m} \leftrightarrow \text{PAS sampling depth } \mu = 7.6 \mu\text{m}$

Photoacoustic signal observed beyond μ

2. PAS spectrum of multilayered film, various mirror velocity


Multilayered polyimide film
Kapton® 300F929



Estimation of accurate sampling depth $A\mu$ (A: factor)

- 1775 cm^{-1} peak: Observed at 10 kHz, not observed at 20 kHz
 - » $3.8A < 12.5 \mu\text{m} < 4.8A$
- 1080 cm^{-1} peak: Observed at 10 kHz, not observed at 20 kHz
 - » $3.4A < 12.5 \mu\text{m} < 4.8A$
- 605 cm^{-1} peak: Observed at 20 kHz, not observed at 40 kHz
 - » $4.6A < 12.5 \mu\text{m} < 6.4A$

Estimated factor
 $2.6 < A < 2.7$
i.e., $A \approx e$



Conclusion

Photoacoustic signal can be detected beyond PAS sampling depth μ .
Accurate sampling depth is estimated to be $\sim e\mu$, where the generated signal decays to $1/e^2$ upon reaching the sample surface.